

SSD-Bedkit Wired Bed Exit Alarm

The risk of a patient getting out of their bed and leaving themselves susceptible to injury is a problem many caregivers have to deal with on a daily basis. We offer a premium quality bed exit alarm that make this job easier and helps reduce the likelihood a patient will leave their bed and wander unsupervised while at risk of falling. The system can be **connected to the nurse call system** using the provided connection lead or can be **used stand-alone** if preferred.



What does it do?

The bed exit alarm notifies staff when a resident leaves their bed and begins to wander. The alarm can sound locally in the room itself or it can be connected to a nurse call system so that the alarm is transmitted around the building. The system can also be supplied with portable pagers so that the alarm is sent to staff directly while they are on the move.

Features:



Available as either bed or chair sensor mat



Optional 2 second delay to help minimise false alarms



Silent mode available when connected to call systems



Clip-on sensor can raise an alarm if client leaves area



Multiple audio levels to suit every situation



Battery-powered with optional mains PSU available

How does it work?

The system is supplied with an alarm box, a bed sensor mat and a nurse call connection lead. The sensor mat is thicker and softer than standard sensor mats and features a robust connection to help prevent damage if the wire was to get pulled. It comes in its own protective card sleeve to further protect the mat from damage.

The sensor mat is placed on top of the bed, underneath the sheets. The alarm box can then be fixed to the bed frame using the provided hanging brackets or placed on the bedside table.

Note: If the facility has a nurse call system, the alarm box can be connected to this by plugging it into a nearby call point so that calls are sent around the building. In facilities without a nurse call system, the alarm box will sound an alert locally in the room!

Contents:

- 1 x Premium quality bed sensor mat (size: 75cm x 25cm)
- 1 x Control box
- 1 x Nurse call connection lead
- 1 x Set of hanging brackets
- 1 x 9V battery

For more information on these products, please contact us: